

AC-1100 COLPOSCOPE

ALLTION AC-1100 colposcopes have two models now. One is AC-1100 with 6V/20W Halogen Light Source with heat insulating house. Another is AC-1110 with 5.9V/5W LED Cold Light Source giving white color illumination.

This is a high quality and most basic colposcope that offers excellent optics to your practice affordable. This unit is compact, light-weight and with roller base for easily maneuver.

ALLTION colposcopes are examined and calibrated prior to shipping, so they arrive completely assembled and ready to use. All you have to do is plug it in.

FEATURES OF ALLTION AC-1100 & AC-1110 COLPOSCOPE:

1. Precision computer designed total magnification 9x
2. 316mm focal length
3. Field of View: 22mm dia.
4. Depth of Field: > 4mm
5. Light intensity: >10,000 LUX
6. Individually adjustable eyepiece 12.5x (diopters)
7. Right eyepiece built-in two circles for you to measure the size of the problem you see.
8. Built-in rheostat for brightness adjustment
9. User-selectable built-in green filter for enhanced contrast
10. The halogen bulb or LED bulb is easily changeable.
11. Working Height: 838mm to 1194mm
12. Two (Gross/Fine) focus adjustment Handles
13. Two (Gross/Fine) height adjustment Handles
14. Universal Teflon[®] ball-joint gives you 360° easy motion
15. Unit is completely assembled and ready for use once power is on
16. Main Power Supply: 110V/220V



Two Circles Above for Measuring Size of Tissue



CE No. 075S019-JT-CE-P02V02

FDA 510K Number: K061306

ISO9001: 2000 Certificate No. 5300/2002137799

*Items and specification subject to change,
please contact ALLTION for current product information*

ALLTION

(Wuzhou) Co., Ltd.

Microscopes Are Our Business & Service Is Our Pride

Address: Suite B, 11/F., Lijing Tower, Cathay Plaza, No. 18, Zhong Shan Road, Wuzhou, Guangxi, China 54300

Tel: +86-774-2836101 Fax: +86-774-2836192

E-mail: sales@alltion.com or sales@alltion-microscope.com

Web site: <http://www.alltion.com> or <http://www.alltion-microscope.com>