

AC-1300 COLPOSCOPE



ALLTION AC-1300 colposcopes have two models, AC-1300 with 6V/20W Halogen Light Source with heat insulating house and AC-1310 with 5.9V/5W LED Cold Light Source giving white color illumination.

This 3-steps magnification colposcope is a high quality basic colposcope that offers excellent optics to your practice affordable. This unit is compact, light-weight and with roller base for easily maneuver.

ALLTION colposcopes are examined and calibrated prior to shipping, so the arrives are completely assembled and ready to use. All you have to do is plug it in.

FEATURES OF ALLTION AC-1300 & AC-1310 COLPOSCOPE:


1. Total magnification 3.75x, 7x, 15x
2. 300mm focal length
3. Field of View: 30.93mm, 15.50mm, 7.74mm dia.
4. Depth of Field: 18mm, 4mm, 2mm
5. Light intensity: >10,000 LUX
6. Individually adjustable eyepiece 12.5x (diopters)
7. Right eyepiece with built-in two circles for you to measure the size of the problem you see.
8. Built-in rheostat for brightness adjustment
9. User-selectable built-in green filter for enhanced contrast
10. The halogen bulb or LED bulb is easily changeable.
11. Working Height: 838mm to 1194mm
12. Two (Gross/Fine) focus adjustment Handles
13. Two (Gross/Fine) height adjustment Handles
14. Universal Teflon® ball-joint gives you 360° easy motion
15. Unit is completely assembled and ready for use once power is on
16. Main Power Supply: 110V/220V



CE No.075S019-IT-CE-P02V02

FDA 510K Number: K061306

ISO9001: 2000 Certificate No. 5300/2002137799



ALLTION