

# AC-3000 SERIES COLPOSCOPES




German Optical  
Glass Inside

ALLTION is proud to offer the AC-3000 Series Colposcopes. The lenses inside these colposcopes are made by Schott® Optical Glass imported from Germany, which offers excellent optic feature, such as high resolution, long depth of field, etc.

## FEATURES OF MODEL AC-3300 COLPOSCOPE

- Smooth three-step turret magnification changer.
- 300mm working distance objective lens offering 2.7x, 6.7x, and 16.7x magnifications. Optional magnifications/objective lens combinations are available.
- Widefield 12.5x eyepieces with diopter locks and fold-down rubber eyecups for eyeglass wearers.
- Straight binocular head F=160.
- Light Intensity on object plane: 38,000.00 Lux. Min.
- Fan cooled, \*\*dual-port variable fiber optic coaxial illumination with 15V/150W lamp  
*\*\*should a lamp fail during an examination, the dual-port feature allows a quick change from one lamp to another resulting in an uninterrupted procedure.*
- Maximum height from floor to objective lens: 1200mm
- Minimum height from floor to objective lens: 1050mm
- Green filter in selector knob.
- Small diameter base on roller castors with adjustable vertical post height.
- Spare lamp, dust cover and instruction manual.
- Video accessories available (see accessory page).



CE

CE No. 076S050-IT-CE-P02V02  
ISO 9001:2000/ISO 13485:1996 No.2003-OSL-AQ-7346 issued by DNV  
ISO 14001:2004 NO. 1305-2004-AE-RGC-RvA issued by DNV  
OHSAS 18001:2004 No. 1306-2004-HSO-RGC-DNV issued by DNV

ALLTION